


**WHAT CAN WE
LEARN ABOUT
DIFFUSION BY
NMR ?**



Content


- 1) Diffusion in NMR
- 2) Convection
- 3) The pulse sequences
- 4) Setup and analysis
- 5) Mutual diffusion and diffusion of pure compounds




What is diffusion?

Ecuación Stokes-Einstein



$$D = \frac{k_b T}{6\pi\eta r_s}$$




Self-diffusion




Pulse field gradients





Diffusion measurement: gradient echo






Convection





Diffusion vs convection




- Convection can be compensated
- Diffusion cannot

Diffusion vs convection





Convection compensated spin echo




- Convection is refocused
- Diffusion isn't

Measurement of convection



Diffusion + convection




$$I_G \approx I_{G=0} \cos(\gamma \delta G \Delta v) \exp\left(-(\gamma \delta G)^2 D \left(\Delta - \frac{\delta}{3}\right)\right)$$

convection diffusion

Convection is almost always present

Convection in liquid-state NMR: expect the unexpected

T. M. Barbosa,^{a,b} R. Rittner,^a C. F. Tormena,^a G. A. Morris,^b and M. Nilsson^b




5mm

3mm

- Convection is almost often present
- it depends on T, solvent, probe, geometry of the tube






PFG and Stimulated Echo



- Diffusion may be slow compared to T_2
- Stimulated Echo protects the magnetization in z
- STE drop half of magnetization but is still more intense

BPPLED and Oneshot





- Bipolar Pulse Pair Longitudinal Eddy Current Delay
 - **reduce eddy current BPP**
 - long phase cycling
- Oneshot
 - **short phase cycling** (one scan)
 - small tilt angle
 - lock refocused by small gradients



First detect convection

Una forma rápida de rastrear la presencia de convección es usar diferentes valores del diffusion delay para una misma temperatura de trabajo.







Relación entre el
Dapp y el Big delta


$$D_{\text{app}} = D + \frac{\Delta v^2}{2}$$

- Convection must be detected first since calibration is performed using known solvents

The calibrate your gradients




Recalculate the gradients



$$S(q) = S(0) \exp \left[-D\gamma^2 g^2 \delta^2 (\Delta - \delta/3 - \tau/2) \right]$$

$$g = \sqrt{-\frac{1}{D\gamma^2 \delta^2 (\Delta - \delta/3 - \tau/2)} \ln \left[\frac{s(q)}{s(0)} \right]}$$



$$g = \sqrt{\frac{1}{D\gamma^2 \delta^2 (\Delta - \delta/3 - \tau/2)} \ln \left[\frac{s(0)}{s(q)} \right]}$$

Minimize convection

Rayleigh Number


1) Tube diameters :

$$R = \frac{g\alpha}{\kappa\nu} r^4 T^*$$

2) Solvent viscosity:

Solvent	Viscosity (mPa.s)
Acetone	0.31
Acetonitrile	0.35
Benzene	0.60
Chloroform	0.54
Dichloromethane	0.41
Dimethylsulfoxide	1.99
Methanol	0.54
Tetrahydrofuran	0.46
Toluene	0.56
Water	0.89

3) Sample spinning :



- Sample spinning rate should be synchronized with Δ

Summary

3mm

S. Spinning

BPPLED

One shot




What NMR can tell us about diffusion?


No force applied



self-diffusion




Intra-diffusion



3 mm
S. Spinning
One shot
BPPLED

Force is applied



Mutual diffusion



Fick vs Maxwell - Stefan




Adolf Eugen Fick



Josef Stefan




James Maxwell



$$\mathfrak{D} = \frac{D}{\Gamma}$$

Thermodynamic
correction factor


Maxwell-Stefan



thermodynamic interaction
vs.
molecular friction

D can be measured, while \mathfrak{D} can be calculated

The route to Fick diffusivities



Darken model


$$\mathcal{D}_{12} = x_2 \hat{\mathcal{D}}_1 + x_1 \hat{\mathcal{D}}_2$$

Vignes model



$$\mathcal{D}_{12} = (\mathcal{D}_{12}^{x_1 \rightarrow 1})^{x_1} (\mathcal{D}_{12}^{x_2 \rightarrow 1})^{x_2}$$




require diffusivities
at infinite diffusion
(large error in predicted ones)




MD simulation requires self-diffusion



Molecular Dynamics



Self diffusion coefficient





UNIVERSITÄT PADERBORN
Die Universität der Informationsgesellschaft

Awful results!


- Convection
- Eddy current
- Radiation damping
- Relaxation times
- lock?

Radiation Damping (RD) / high concentrated samples


- dilution
- sample volume




Q-switch



Small tilt angle



Effects of RD




<http://anorganik.uni-tuebingen.de/klaus/nmr/index.php?p=damping/damping>

Effects of RD on D

High- Q




Q -switched




100 Hz

(A)



(B)





Concentric tubes

Anal. Chem. 2002, 74, 191–198



Quantitative NMR Spectroscopy Using Coaxial Inserts Containing a Reference Standard: Purity Determinations for Military Nerve Agents

Terry J . Henderson*

Edgewood Chemical–Biological Forensic Analytical Center, Battelle Memorial Institute Edgewood Operations,
2012 Tollgate Road, Suite 206, Bel Air, Maryland 21015



Accurate determination of diffusion of pure water




- Holz, 2000
- NMR_coaxial_tube

System: H ₂ O 99.9% / DMSO-d ₆		
Temperature (K)	D (m ² /s)	% Error
298	2.257	1.83
308	2.900	-0.17
318	3.614	-0.36
328	4.297	0.83


One shot 15° S. Spinning T.C. 3mm

- Holz measured at low field (>RD) with special probe (>convection)


Accurate determination of diffusion of cyclohexane



Accurate determination of diffusion of lutidine



$$\frac{I(g)}{I(0)} = \exp\{-\gamma^2 g^2 \delta^2 [\Delta + \delta(\alpha^2 - 2)/6 + \tau(\alpha^2 - 1)/2]\}$$



□ (ppm)

One shot




15°

S. Spinning

T.C.


3mm

DOSY representation of D of lutidine




General NMR Analysis Toolbox


What NMR can tell us about diffusion?




Force is applied




Experimental set-up



Results - time evolution of gradients



Results



Darken-Schneider



$$D_{12} = (x_2 D_1^* + x_1 D_2^*) \left[1 + \frac{d \ln \Gamma_1}{d \ln x_1} \right]$$

D'Agostino (2012)



$$D_{12} = (x_2 D_1^* + 2x_1 D_2^*) \left[1 + \frac{d \ln \Gamma_1}{d \ln x_1} \right]^\alpha$$

dimerization

Clustering factor



What NMR can tell us about diffusion?



Acknowledgements



Christian
Pantoja



Lorraine
Guascar



UNIVERSITÄT PADERBORN
Die Universität der Informationsgesellschaft

REVIEW

[View Article Online](#)
[View Journal](#)



Cite this: DOI: 10.1039/c7np00064b

The value of universally available raw NMR data for transparency, reproducibility, and integrity in natural product research[†]

James B. McAlpine,^{1,2*} Shao-Nong Chen,^{1,3} Andrei Kutateladze,^{1,4} John B. MacMillan,^{1,5} Giovanni Appendino,^{1,6} Andersson Barison,^{1,7} Mehdi A. Beniddir,^{1,8} Maike W. Biavatti,^{1,9} Stefan Bluml,^{1,10} Asmaa Boufridi,^{1,11} Mark S. Butler,^{1,12} Robert J. Capon,^{1,13} Young H. Choi,^{1,14} David Coppage,^{1,15} Phillip Crews,^{1,16} Michael T. Crimmins,^{1,17} Marie Csete,^{1,18} Pradeep Dewapriya,^{1,19} Joseph M. Egan,^{1,20} Mary J. Garson,^{1,21} Grégory Genta-Jouve,^{1,22} William H. Gerwick,^{1,23} Harald Gross,^{1,24} Mary Kay Harper,^{1,25} Precilia Hermanto,^{1,26} James M. Hook,^{1,27} Luke Hunter,^{1,28} Damien Jeannerat,^{1,29} Nai-Yun Ji,^{1,30} Tyler A. Johnson,^{1,31} David G. I. Kingston,^{1,32} Hiroyuki Koshino,^{1,33} Hsiau-Wei Lee,^{1,34} Guy Lewin,^{1,35} Jie Li,^{1,36} Roger G. Linington,^{1,37} MiaoMiao Liu,^{1,38} Kerry L. McPhail,^{1,39} Tadeusz F. Molinski,^{1,40} Bradley S. Moore,^{1,41} Joo-Won Nam,^{1,42} Ram P. Neupane,^{1,43} Matthias Niemitz,^{1,44} Jean-Marc Nuzillard,^{1,45} Nicholas H. Oberlies,^{1,46} Fernanda M. M. Ocampos,^{1,47} Guohui Pan,^{1,48} Ronald J. Quinn,^{1,49} D. Sai Reddy,^{1,50} Jean-Hugues Renault,^{1,51} José Rivera-Chávez,^{1,52} Wolfgang Robien,^{1,53} Carla M. Saunders,^{1,54} Thomas J. Schmidt,^{1,55} Christoph Seger,^{1,56} Ben Shen,^{1,57} Christoph Steinbeck,^{1,58} Hermann Stuppner,^{1,59} Sonja Sturm,^{1,60} Orazio Tagliaferla-Scafati,^{1,61} Dean J. Tantillo,^{1,62} Robert Verpoorte,^{1,63} Bin-Gui Wang,^{1,64} Craig M. Williams,^{1,65} Philip G. Williams,^{1,66} Julien Wist,^{1,67} Jian-Min Yue,^{1,68} Chen Zhang,^{1,69} Zhengren Xu,^{1,70} Charlotte Simmler,^{1,71} David C. Lankin,^{1,72} Jonathan Bisson^{1,73} and Guido F. Pauli^{1,74*}

FEATURE ARTICLE RESEARCH ARTICLE

NMReDATA, a standard to report the NMR assignment and parameters of organic compounds

Marion Pupier¹ | Jean-Marc Nuzillard² | Julien Wist³ | Nils E. Schlörer⁴ | Stefan Kuhn⁴ | Mate Erdélyi⁵ | Christoph Steinbeck⁶ | Antony J. Williams⁷ | Craig Butts⁸ | Tim D.W. Claridge⁹ | Bozhana Mikhova¹⁰ | Wolfgang Robien¹¹ | Hesam Dashti¹² | Hamid R. Eghbalian¹² | Christophe Farès¹³ | Christian Adam¹⁴ | Pavel Kessler¹⁵ | Fabrice Moriaud¹⁶ | Mikhail Elyashberg¹⁷ | Dimitris Argyropoulos¹⁸ | Manuel Pérez¹⁹ | Patrick Giraudeau^{20,21} | Roberto R. Gil²² | Paul Trevorrow²³ | Damien Jeannerat¹



SPECIAL ISSUE RESEARCH ARTICLE | Open Access

The C6H6 NMR repository: An integral solution to control the flow of your data from the magnet to the public

Luc Patiny, Michaël Zasso, Daniel Kostro, Andrés Bernal, Andrés M. Castillo, Alejandro Bolaños, Miguel A. Asencio, Norman Pellet, Matthew Todd, Nils Schloerer, Stefan Kuhn, Elaine Holmes, Sacha Javor, Julien Wist✉, ... See fewer authors ^

First published: 05 October 2017 | <https://doi.org/10.1002/mrc.4669>

Welcome to www.cheminfo.org

From this webpage you will be able to try many applications using the open source project [Visualizer](#).

You may be interested in trying those fancy tools:

- . [Wikipedia chemical structure explorer](#)
- . [PDB explorer](#)
- . [Virtual combinatorial library](#)
- . [Electronic laboratory notebook \(ELN\) : www.c6h6.org](#)

Want to create tools by yourself ?

- . Just LOGIN (using a Gmail account) and you will be able to edit and save the views from the "View Manager" button.
- . Check simple views in "Tutorial -> Introduction"

You may also jump to specific selection of tools:

Mass spectrometry

NMR spectroscopy

Student / teacher in:

- General chemistry
- Organic chemistry
- Biooriented chemistry
- Structural analysis (NMR, IR, ms)
- Cheminformatics

<http://www.nmrdb.org>

<http://www.c6h6.org>



ÉCOLE POLYTECHNIQUE
FÉDÉRALE DE LAUSANNE



Universidad
del Valle

